

**Safeguarding Hygiene
Where it Matters Most**

HyServe



Lumitester™ Smart

Hygiene monitoring
in 10 seconds



**HYGIENE
MONITORING IN
THREE EASY STEPS**



**VISUALIZE THE
DATA VIA APP**



**ACCESS FROM
ANYWHERE WITH
CLOUD STORAGE**

+ LucIPac™
With A3 Detection

ATP + ADP + AMP

More Reliable and Sensitive than ATP Tests
A3 Finds What Others Miss



MAKE HYGIENE MONITORING A HABIT FOR EVERYONE

With the App, anyone can check cleanliness easily



Measure

Easy-to-use, results within 10 seconds



Swab

Swab the sample on a 10 x 10cm area, and insert the swab into main body



Measure

Shake well and insert into the Lumitester



Analyze

Measurements (RLU) displayed within 10 seconds



Understand

App to automatically analyse data

Recorded data can be displayed in time-series by registering with the Lumitester Smart App. Inspection pass rates are automatically graphed and data can be visualized. Employees will be more conscious of hygiene, and will maintain a high standard of environmental sanitation, which will help build greater trust in your business.



Display trend analysis for each inspection



Visualize overall inspection scores through graphs



Connect

Cloud data management to combine data from multiple facilities

Data will be stored in the cloud and can be accessed from anywhere in the world. This will allow fast detection and response to any problems that may arise.



Lumitester™ Smart

Product code: 1902653

Measurement time	10 seconds
Detection	ATP + ADP + AMP
Detection method	Luciferase-Enzyme Assay
Unit of Measurement	Relative Light Units (RLU)
Included Accessories	2 alkaline batteries (AA), Cleaning brush, USB cable, Strap, Quick Manual

- Use for hygiene monitoring and other applications where detection of organic residues needs to be performed.

- Lucipac A3 and Lumitester Smart should not be used for counting general living bacteria or detecting specific pathogens.

LuciPac™ A3 Surface

Product code: 1702671 100 swabs/kit

LuciPac™ A3 Water

Product code: 1702672 100 swabs/kit

Long Term Storage	2-8°C (do not freeze)
Short Term Storage	14 days at 25°C (unopened)
Shipping Condition	at ambient temperature for up to 14 days
Expiry	15 months after manufacturing date

distribuito da:



www.bioapp.it
via L.Einaudi 127 - 45100 Rovigo - Italy
0425 1547882 - info@bioapp.it

Use LuciPac A3 for Lumitester Smart, PD-30 or PD-20. Do not use it for other models.