

CryoDryer 5

Ice condenser capacity: 5Kg
 Ice condenser performance: 4Kg/24hr
 Ice condenser temperature: -50°C
 Chamber volume: 65L
 Chamber dimension: 370 x 600 mm
 Shelf dimensions: 535 x 165 mm (B x T)
 Number of shelves: 6 + 1
 Shelf area: 0,529 m²
 Shelf spacing: 19 mm - 217 mm
 Shelf temperature: +40°C (up to 60°C)
 Material: st. steel 304 or 316
 Shelf cooling system: No shelves cooling

Accessories:

- Electromagnetic pressure control and stop valve
- Manifold horizontal with rubber valves for connection of 4 round flasks



CryoDryer 10

Ice condenser capacity: 10Kg
 Ice condenser performance: 8Kg/24hr
 Ice condenser temperature: -50°C
 Chamber volume: 100L
 Chamber dimension: 367 x 508 x 600 mm
 Shelf dimensions: 535 x 330 (B x T)
 Number of shelves: 6 + 1
 Shelf area: 1,058 m²
 Shelf spacing: 19 mm - 217 mm
 Shelf temperature: +40°C (up to 60°C)
 Material: st. steel 304 or 316
 Shelf cooling system: No shelves cooling

Accessories:

- Electromagnetic pressure control and stop valve
- Manifold horizontal with rubber valves for connection of 4 round flasks



Subject to change without notice.
 All information in this brochure is purely informative and non-binding.
 Our quotations are authoritative with regard to orders.



distribuito da:

bioapp
 Rovigo
 YOU THINK, WE DO!

www.bioapp.it

via L.Einaudi 127 - 45100 Rovigo - Italy

0425 1547882 - info@bioapp.it