




CryoDryer 100

Ice condenser capacity: 100Kg
 Ice condenser performance: 80Kg / 24hr
 Ice condenser temperature: -50°C / -80°C
 Chamber volume: 3000 L
 Shelf dimensions: 800 x 2000 mm (B x T)
 Number of shelves: 7 + 1
 Shelf area: 10 m²
 Shelf spacing: 40 mm
 Shelf temperature: -35°C to +40°C for ice condenser at -50°C
 Shelf temperature accuracy: max. ±2 °C
 Ice condenser cooling system: Air cooled (optionally water cooled)
 Shelf cooling system: (see ice condenser cooling system)
 Chamber dimensions: 1300 x 2400 mm (diam x D)

CryoDryer 300

Ice condenser capacity: 300Kg
 Ice condenser performance: 280Kg / 24hr
 Ice condenser temperature: -50°C / -80°C
 Chamber volume: 7500 L
 Shelf dimensions: 1200 x 2000 mm (B x T)
 Number of shelves: 12 + 1
 Shelf area: 29 m²
 Shelf spacing: 40 mm
 Shelf temperature: -35°C to +40°C for ice condenser at -50°C
 Shelf temperature accuracy: max. ±2°C
 Ice condenser cooling system: Air cooled (optionally water cooled)
 Shelf cooling system: (see ice condenser cooling system)
 Chamber dimensions: 2000 x 3000 mm (diam x D)


CryoDryer 600

Ice condenser capacity: 600Kg
 Ice condenser performance: 580Kg / 24hr
 Ice condenser temperature: -50°C / -80°C
 Chamber volume: 15000 L
 Shelf dimensions: 1200 x 4000 mm (B x T)
 Number of shelves: 12 + 1
 Shelf area: 58 m²
 Shelf spacing: 40 mm
 Shelf temperature: -35°C to +40°C for ice condenser at -50°C
 Shelf temperature accuracy: max. ±2°C
 Ice condenser cooling system: Air cooled (optionally water cooled)
 Shelf cooling system: (see ice condenser cooling system)
 Chamber dimensions: 2000 x 6000 mm (diam x D)



CryoDryer 20

Ice condenser capacity: 20Kg
Ice condenser performance: 18Kg / 24hr
Ice condenser temperature: -50°C
Chamber volume: 300L
Shelf dimensions: 500 x 700 mm (B x T)
Number of shelves: 6 + 1
Shelf area: 1,95 m²
Shelf spacing: 40 mm
Shelf temperature: -35°C to +40°C
Shelf temperature accuracy: max. ±2°C
Ice condenser cooling system: Air cooled (optionally water cooled)
Shelf cooling system: (see ice condenser cooling system)
Unit dimensions: 1878 x 1082 x 1060 mm (H x W x D)



Subject to change without notice.
All information in this brochure is purely informative and non-binding.
Our quotations are authoritative with regard to orders.

distribuito da:

bioapp
Rovigo
YOU THINK, WE DO!

www.bioapp.it

via L.Einaudi 127 - 45100 Rovigo - Italy

0425 1547882 - info@bioapp.it